




Inner Inside
Make 2 - .25 thick
Cut slots .125 deep




Inner Outside
Make 2 - .25 thick
Cut slots .125 deep




Center Rib
Make 2 - .125 thick




Middle Rib
Make 4 - .125 thick




Edge Rib
Make 4 - .125 thick




Ring
Make 2 - .125 thick



Ring End
Make 2 - .125 thick



End Curve
Make 2 - .125 thick
Chamfer End at 30 degrees



Shelf
Make 2 - .125 thick
Cut slots .0625 deep




Bottom Skin
Make 1 - .040 thick



Side Skin
Make 2 - .040 thick



Shelf Skin
Make 2 - .040 thick



Curve Skin
Make 2 - .040 thick

Skin parts are oversize, to be trimmed to fit after being glued in place with the exception of the side skin which is cut to exact width. See assembly instructions for more information.

Eggcrate Center Ankle Page 1/2

Copyright Media-Conversions - Released for non-commercial use ONLY!

DRAWN	DATE				
CHECKED	6/12/2012				
QA					
MFG					
APPROVED		SIZE	FSCM NO.	DWG NO.	REV
		C		Center Ankle	
		SCALE			SHEET